

inspired™



Scan for more information.

# inspired O2FLO™

High Flow Oxygen Therapy is a simple-to-use system that delivers warm, moist gas at high flow rates that generates positive airway pressure. When used at flow rates of 1- 2 L/kg/min it acts as a bridge between Low Flow Oxygen Therapies and Continuous Positive Airways Pressure (CPAP), reducing the need for intubation.



VHC60 Autofeed Humidification Chamber



Nasal Cannula (Adult)

Model	Flow Rate (LPM)
VANC-01 (S)	10-50
VANC-02 (M)	10-60
VANC-03 (L)	10-60



510-049/050  
Integrated Heated  
Breathing Circuit  
(Built-in Temperature Sensor)



Nasal Cannula (Infant)

Model	Flow Rate (LPM)
VINC-01 (S)	2~8
VINC-02 (M)	2~8
VINC-03 (L)	2~10
VINC-03 (XL)	2~10





## O2FLO High Flow Heated Respiratory Humidifier (VUN-001)

Flow Rate Range	10~60 LPM (Adult Mode) 2~25 LPM (Paediatric Mode)
Humidification Output	Complies with ISO 80601-2-74 standard
Dew Point Temperature Range	31- 37 °C, ±2 °C (Adult Mode) 31- 34 °C, ±2 °C (Paediatric Mode)
Oxygen Concentration Range	21%-95%, Accuracy ±5%
Sound Pressure Level	≤50 dB(A)
Power Input	100-240V~, 50/60Hz
Net Weight	2.8±0.2kg
Size	243mm x 220mm x 170mm
Standard	IEC/EN 60601-1, IEC/EN 60601-1-2, ISO 80601-2-74

# ACCESSORIES LIST

Adult 51006729  
Pediatric 51006728




510-090/091  
Heated Breathing Circuit  
10 pcs/box

51006737



VUN-200  
Temperature Probe  
(Used for 510-090/091)

51006736



Heater Connection Cable  
(Used for 510-090/091)

Adult 51006765  
Pediatric 51006766



510-049/050  
Integrated Heated  
Breathing Circuit  
(Built-in Temperature Sensor)  
10 pcs/box

51006767



VUN-300  
Connection Cable  
(Used for 510-049/050)

VHC60



Autofeed  
Humidification Chamber  
10 pcs/box

VANC-01(S) VANC-02 (M)  
VANC-03 (L)



Nasal Cannula (Adult)  
20 pcs/box

VINC-01 (S) VINC-02 (M)  
VINC-03 (L) VINC-04 (XL)



Nasal Cannula (Infant)  
200 pcs/inner box

51006764




Filter  
2000 pcs/box

51006741



Oxygen Tube  
100 pcs/box



**Flow Meter**

We offer a wide range of direct-probes to meet different international standards. Please contact sales representative for model specific to your requirements.



**Trolley**

We offer several options for Trolleys. Please contact sales representative for more details.

## Headquarters

Flat/RM B2, 7/F, Hang Fung Industrial Building, Phase 2,  
2G Hok Yuen Street, Hung Hom, Kowloon, Hong Kong, China

+852 2365 5688 +852 2765 8428

info@inspired-medical.com

## Sales Office-China

1st Floor, Building 10, Guantai Biotech Technical Cooperation  
Center, Songshan Lake Science and Technology Industrial  
Park, Dongguan Guangdong 523808, China

+86 769 2166 7698

## Sales Office-Japan

302, NAKAGAWA Bldg. 3-31-1 Yushima,  
Bunkyo-ku, Tokyo 113-0034, Japan

+81 3 5817 8471

## Service support:

+852 2186 1010

service@inspired-medical.com



www.inspired-medical.com